


Healthy peristomal skin: a foundation for successful voice and pulmonary rehabilitation



Abstract

Heat and Moisture Exchangers (HMEs) and their attachment devices are important medical devices for pulmonary and voice rehabilitation after total laryngectomy. Their positive impact is evident, particularly when adhering to HME use. Adherence to HME use and efficiency of tracheoesophageal (TE) voice production can be challenged by Medical Adhesive-Related Skin Injuries (MARS!) and Peristomal Skin Complications (PSC). Understanding of the skin and the factors that influence peristomal skin and skin in general after total laryngectomy helps identify patients that may be at risk for development of MARS!. Availability of new peristomal adhesives based on modern adhesive materials and technology, careful adhesive selection and appropriate application and removal techniques can help prevent MARS! and maintain the integrity of the peristomal skin. Healthy peristomal skin forms the foundation for adherent HME use and an airtight seal for efficient TE speech.

Table of Contents

1	Background
	1.1 Attachment of Heat and Moisture Exchangers
	1.2 The challenge with peristomal medical adhesives
	1.3 Medical Adhesive-Related Skin Injuries
2	Research for a deeper insight
	2.1 A wide variety of stoma configurations
	2.2 Importance of peristomal adhesives for adherence to HMEs
	2.3 Relationship between adhesive problems and skin irritation
	2.4 Factors influencing adhesive problems
	2.5 Factors influencing peristomal skin irritation
	2.6 Peristomal adhesives and voice rehabilitation
3	Medical adhesives for TL patients – technology
	3.1 General composition of medical adhesives
	3.2 Material characteristics
	3.3 Design characteristics
4	The normal skin
	4.1 The structure of the skin
	4.1.1 The epidermis
	4.1.2 The dermis
	4.1.3 The subcutaneous tissue
5	Medical adhesive-related skin injuries (MARSI)
	5.1 Categories of MARSI
	5.2 Risk factors, causes, and prevention
	5.2.1 Intrinsic and extrinsic risk factors
	5.2.2 Preventable causes of MARSI
	5.2.3 Medical adhesive-related challenges specific for TL patients
6	Provox® Life™ Adhesives
	6.1 Material characteristics and technical details
	6.2 Shape
	6.3 Connector
7	Clinical evidence
	7.1 Post-market User Experience Program
	7.2 Post-market clinical follow-up studies
8	Prevention of MARSI in TL patients
9	Conclusions
10	Q&A with the expert

List of abbreviations

ASL	Airway Surface Liquid
CASA-Q	Cough And Sputum Assessment Questionnaire
HCP	Healthcare Professional
HME	Heat and Moisture Exchanger
MARSI	Medical Adhesive Related Skin Injuries
MTR	Moisture Transport Rate
PSC	Peristomal Skin Complications
QoL	Quality of Life
SCM	Sternocleidomastoid muscle
TE	Tracheoesophageal
TL	Total Laryngectomy
VP	Voice Prosthesis

1. Background

1.1 Attachment of Heat and Moisture Exchangers

Heat and Moisture Exchangers (HMEs) are essential for optimal pulmonary rehabilitation after total laryngectomy (TL) and their benefits have been extensively documented⁽¹⁻⁵⁾. They are usually attached in front of the tracheostoma with a peristomal adhesive, although some patients require a laryngectomy tube or stoma button to prevent narrowing of the tracheostoma. The use of these intratracheal devices can pose problems such as stoma irritation leading to further narrowing of the stoma, crust formation and more serious consequences such as mucosal infection, difficulty clearing secretions and obstruction^(6, 7). Furthermore, air escape around the tube complicates voice rehabilitation and the use of heat and moisture exchangers⁽⁷⁾. The COVID-19 pandemic further emphasized the importance of a complete seal provided by peristomal adhesives; when properly applied they prevent airflow outside of the HME and reduce mucus contamination of clothes and other physical barriers⁽⁸⁾. Furthermore, it was demonstrated that TL patients using laryngectomy tubes pose a higher risk for spreading viral aerosols due to increased velocity of exhaled airflow caused by the presence of the tube as a narrowing factor; it is therefore recommended to use an HME with viral filter and preferably to change laryngectomy tubes to peristomal adhesives when possible⁽⁹⁾. Using HMEs together with an adhesive provides an airtight seal which makes the conditioning of air more reliable and facilitates speech in patients utilizing voice prostheses (VPs)⁽³⁾.

Using HMEs together with an adhesive provides an airtight seal which makes the conditioning of air more reliable and facilitates speech in patients utilizing voice prostheses

1.2 The challenge with peristomal medical adhesives

Patient adherence to HME use has been identified as an issue and a large variation (35%-73%) has been reported⁽¹⁰⁻¹³⁾. Intolerance to the breathing resistance is one reason for reduced adherence to HME use and this has been addressed with the development of Provox® Life™ HMEs and is covered in detail in a previous Whitepaper by Atos Medical⁽¹⁴⁾. Another common reason for reduced adherence to HME use is problems associated with HME peristomal adhesive baseplates and other attachments. These issues include skin irritation, poor attachment, as well as poor fit to stoma shape, particularly for the deeper stomas^(1, 13, 15).

A common reason for reduced adherence to HME use is problems associated with HME peristomal adhesive baseplates and other attachments

Although deep stomas may be prevented to some extent by applying surgical techniques such as cutting of the heads of the Sternocleidomastoid (SCM) muscles, horizontal transection and suturing techniques preventing exposure of cartilage⁽¹⁶⁾, prioritization of oncological outcomes and the need for flap reconstructions can lead to difficult stomas. Problems with peristomal adhesives are exacerbated in patients using automatic handsfree speaking valves; regular use is only achieved by about 20%-40% of patients, mainly due to problems achieving durable and airtight fixation with the adhesive due to the phonation pressures exerted onto the adhesive during handsfree speech⁽¹⁷⁻¹⁹⁾. Dermatological problems from adhesives have been reported for 20-60% of total laryngectomy patients^(4,13,20) and lead to poor adherence to HME use, especially at night when they frequently take 'skin rest'⁽²¹⁾. The development of more skin-friendly and better fitting alternatives is therefore a promising measure to improve adherence to HME use.

1.3 Medical Adhesive-Related Injuries

Laryngectomy patients are susceptible to Peristomal Skin Complications (PSC) and Medical Adhesive-Related Injuries (MARSI), just like patients with ostomies or urostomies. A consensus summit on MARSI held in December 2012 with a panel of 23 key opinion leaders concluded that skin care, including its protection against MARSI, is a basic requirement for patient care⁽²²⁾. Although tracheostomas are not mentioned specifically in that paper, the specialty practice areas represented by the members of the panel include closely related areas such as dermatology, oncology, perioperative, plastic surgery, wound, ostomy and continence. Therefore, this consensus work is to a large extent also applicable to the use of peristomal adhesives post-laryngectomy. The expert opinions presented in the paper suggest that appropriate selection of adhesive products and the use of proper application and removal techniques can help minimize MARSI.

This Whitepaper aims to raise awareness for PSC and MARSI in TL patients with healthcare professionals (HCP) providing services to this patient population and to provide knowledge and tools necessary for preventing and managing MARSI and PSC in TL patients.

2. Research for a deeper insight

Published research indicates that issues related to HMEs and attachments limit the possibilities for adherence to HME use among total laryngectomy patients, which could have a negative impact on their Quality of Life (QoL) ^(1, 2, 3, 15). For a deeper insight into the lives and needs of patients in a contemporary context, and a better understanding of issues related to HME and attachment products, Atos Medical initiated an extensive collaborative research project in 2015, which was completed in 2019. Parts of the results were published in Leemans et al. ⁽²³⁾. The research was complemented by a large market research project into HME attachments conducted by Atos Medical in 2022.

RESEARCH METHODS

Stage 1.

Anthropological studies in collaboration with ReD Associates (www.redassociates.com, Copenhagen, Denmark):

- Extensive first study of 26 total laryngectomy patients and >50 health care professionals (HCPs) in four countries (US, Spain, Germany and Sweden) by way of interviews and observations. Data included 100+ hours of video recordings, 1000+ pages of notes and 1000+ photographs.
- In-depth second study focused on pulmonary health and HMEs including 18 patients and 18 HCPs by way of interviews and observations.

Stage 2.

Validation of stage 1 findings by a questionnaire-based quantitative survey in collaboration with ReD Associates and the Netherlands Cancer Institute (Amsterdam, the Netherlands):

- Survey questionnaire covered patient demographics, pulmonary health, impact of total laryngectomy on QoL and use of medical products related to total laryngectomy.
- 1734 total laryngectomy patients in 9 countries (UK, US, Germany, Sweden, France, the Netherlands, Brazil, Italy and Spain) completed the questionnaire (participation rate 21%).
- The respondents' group was demographically representative of the total laryngectomy population and included patients at different postoperative stages using different medical products.

Stage 3.

Market research for better understanding of issues with HME attachments

- A web-based survey questionnaire was sent via e-mail to 16,517 TL patients who were customers of Atos Medical, had an email address on file, and who had consented to receive marketing communications from Atos Medical.
- 2852 TL patients in 12 countries (US, UK, Germany, Italy, France, Canada, the Netherlands, Australia, Portugal, Denmark, Switzerland, New Zealand) completed the survey (participation rate 17%).
- The respondents' group was demographically representative of the total laryngectomy population and included patients of different ages and at different postoperative stages.
- 256 patients did not use any attachment (adhesive, tube, or button) at the time of the survey and were excluded from further analyses, leaving a total of 2596.
- The majority used an adhesive (n=1882, 72.5%) and 27.5% (n=714) used a tube or button. Thirty-eight percent (n=709) of the adhesive users sometimes used a laryngectomy tube or stoma button in addition to their adhesive.

2.1 A wide variety of stoma configurations

The ethnographic research completed during Stage 1 demonstrated that a TL rarely results in a ‘perfect stoma’ due to tissue scarring, reconstructive surgeries and the surgical priority of cancer removal. Thus, the TL patient population exhibits a wide variety of tracheostoma configurations (see Figure 1). HCP wanted their patients to wear HMEs but recognized there are barriers. Adhesive-related problems experienced were adhesives not sticking to the skin, not properly covering the stoma, not forming a seal, irritating the skin, being uncomfortable to wear for a long period of time, and not having the right shape for the stoma. Patients reported red, sore, and itchy skin from adhesive materials and pressure.




FIGURE 1. Collage of laryngectomy stomas, demonstrating the large variety in skin types and stoma configurations.

2.2 Importance of peristomal adhesives for adherence to HMEs

The quantitative survey conducted in Stage 2 provided in-depth insights into problems specifically reported for HME attachments. The results indicate that problems with adhesives play an important role in reduced adherence to HME use. Half of the patients reported that wearing adhesives caused their stomas to become red, raw, or irritated and 48% kept their stoma uncovered to allow the skin to heal or breathe (Table 1). Parts of this research have been published by Leemans et al. (23).

TABLE 1. Percentages of patients reporting problems with HME attachments

HME attachment	Problem with attachment	% of patients
Adhesives	The skin around my stoma has become red, raw or irritated after wearing adhesives	50%
	I have left my stoma uncovered without an adhesive to let my skin heal or breathe	48%
	I have had difficulties getting adhesives to stay on because of my stoma or neck shape	38%
	I have had difficulties getting adhesives to stay on because of sensitive, irritated or wet skin	33%
	My adhesive has come loose or has leaked air earlier than expected	53%
Buttons, tubes	My adhesive, tube or button has come loose because of coughing	51%
Interface	My HME has popped out because of coughing	45%

2.3 Relationship between adhesive problems and skin irritation

The HME attachments market research conducted in Stage 3 indicated that the majority of the 1882 adhesive users used adhesives without any problems (n=1182, 63%), 36% (n=677) experienced problems (with 64 of them experiencing ‘a lot of problems’) and 1% (n=23) responded ‘I don’t know’. Table 2 outlines the specific problems patients most reported.

TABLE 2. Frequency of patient-reported adhesive problems.

Problem with adhesive*	Number (%) of respondents
It did not stick to my skin for as long as I needed it to	275 (41%)
It didn’t provide an airtight seal when speaking	244 (36%)
My skin could not tolerate it (irritation or reaction)	94 (14%)
It did not provide a comfortable fit to my stoma	63 (9%)

* Question (to patients experiencing adhesive problems): Can you please clarify the problems you experience wearing an adhesive? N = 677 patients. Multiple answers possible. Answer options also included a few other options that are not shown here due to low response rates and/or not medical in nature such as reimbursement restrictions, stoma stenosis and not using an HME at the moment.

A total of 1423 patients used tubes and/or buttons of which 714 used them exclusively and 709 used them in combination with an adhesive. Tubes and buttons were mainly used to keep the stoma stable or open (N= 871, 61%), for ‘skin rest’ (N=237, 17%) and for problems with adhesives (N= 175, 12%). Other reasons were ‘preference’, ‘need to use for other reason’, and ‘don’t know’. Tube/button use was more common in patients shorter than 2 years post-surgery with 24% (n=64) of patients < 2yrs post-surgery using a tube or button, versus 16% (n=422) of those 2-6 years post-surgery and 14% (n=352) of those 6 years or longer.

When specifically asked about the frequency and severity of skin irritation, 66% (N=1882) reported this to be 'never' or 'rarely the case, whilst 18% experienced skin irritation 'sometimes', 9% 'often', 5% 'always' and 3% answered 'don't know'. Of the patients with skin irritation, eighty-seven percent described the irritation as mild or somewhat irritated, 7% as very irritated (red, hot, itchy) and 3% as very irritated (breaks down, oozes fluids).

A significant relationship was found between the frequency of peristomal skin irritation and problems with the adhesive (chi-square test, $p < .00001$, see Table 3).

TABLE 3. Cross tabulation of relationship between peristomal skin irritation and adhesive problems (chi-square 126.99, $p < .00001$).

		Skin irritation			Total
		Frequent (always/often)	Infrequent (sometimes)	No skin irritation (rarely/never)	
Adhesive problems	No problems (without any problems)	102 (38.8 %)	182 (52.1 %)	874 (72.1 %)	1158
	Problems (some/a lot of problems)	161 (61.2 %)	167 (47.9 %)	339 (27.9 %)	667
Total		263 (100 %)	349 (100 %)	1213 (100%)	1825

While the adhesive may be the initial cause of the skin irritation, the irritated skin causes problems with the adherence of the adhesive, causing the adhesive to detach, increasing the frequency of replacement, and further exacerbating the problem (see Figure 2). This strong relationship highlights the importance of preventing and managing skin irritation and has been acknowledged in related clinical fields ⁽²⁴⁻²⁶⁾.




FIGURE 2. Vicious circle of adhesive problems leading to skin irritation and irritated skin causing further problems with the adhesive sticking to the skin.

2.4 Factors influencing adhesive problems

As a next step, relationships were explored to further our learnings about adhesive problems, skin irritation and factors of potential influence such as age, postoperative stage, prior radio and/or chemotherapy, stoma depth, neck distension during TE speech, perspiration levels and communication method. These relationships are shown in Table 4.

TABLE 4. Exploration of relationships (chi-squared tests, p-values are shown in the cells) between adhesive problems and age, postoperative stage, (chemo)radiotherapy, stoma depth, speaking method, neck distension and perspiration represented in the columns. n.s.: not significant.

Adhesive problem	
Age	n.s. (.478)
Postoperative stage	.005
CRT	n.s. (.331)
Stoma depth	< .00001
Speaking method	.003
Neck distension during TE speech	< .00001
Perspiration	n.s. (.056)

As can be seen from Table 4, adhesive problems demonstrated a significant relationship with postoperative stage, stoma depth, speaking method, and neck distension during TE speech. Further details on the significant relationships are provided in Table 5.

TABLE 5. Overview of factors significantly related to problems with the adhesive.

		Problems with Adhesive		Level of significance (p)
		No (%)	Yes (%)	
Time since TL surgery (yrs)	<2	59.7	40.3	$p = .005$
	2 – 6	63.4	36.6	
	>6	68.7	31.3	
Stoma depth	Flat	72.3	27.7	$p < .00001$
	Deep	58.2	41.8	
	Very Deep	47.7	52.3	
Speaking method	Voice prosthesis	61.5	38.5	$p = .003$
	No voice prosthesis	68.8	31.2	
Neck distension	No/minor	69.5	30.5	$p < 0.00001$
	Moderate/Much	53.3	46.5	

The details provided in Table 5 indicate that problems with the adhesive are more frequent in the earlier postoperative stages, in patients with deeper set stomas, in patients speaking with a VP, and in patients with neck distension during TE speech. Understanding these factors can help predict which patients are more likely to experience adhesive problems and identify those that may need extra support. Furthermore, surgical techniques aimed at creating flatter stomas should be considered.

2.5 Factors influencing peristomal skin irritation.

The same factors were explored to expand our understanding of peristomal skin irritation and

potential influencing factors such as age, postoperative stage, prior radio and/or chemotherapy, stoma depth, speaking method, neck distension during TE speech and perspiration levels. These relationships are shown in Table 6.

TABLE 6. Exploration of relationships (chi-squared tests, p-values are shown in the cells) between peristomal skin irritation and age, postoperative stage, (chemo)radiotherapy, stoma depth, speaking method, neck distension and perspiration represented in the columns. n.s.: not significant.

Frequency Skin irritation	
Age	.001
Postop	.0001
CRT	n.s. (.960)
Stoma depth	.001
Speaking method	n.s. (.248)
Neck distension during TE speech	.00003
Perspiration	n.s. (0.056)

The frequency of skin irritation showed a significant relationship with age, postoperative stage, stoma depth, neck distension during TE speech and perspiration. Details on the significant relationships are provided in Table 7.

TABLE 7. Overview of factors significantly related to skin irritation.

		Frequency of skin irritation			Level of significance (p)
		Never/Rarely (%)	Sometimes (%)	Often/Always (%)	
Age (yrs)	<39	44.8	22.2	33.3	p = .001
	40-49	51.9	24.1	24.1	
	50-59	60.0	20.7	19.3	
	60+	68.5	18.6	12.9	
Time since TI surgery (yrs)	<2	62.5	19.4	18.2	p = .0001
	2 – 6	64.8	20.2	15.0	
	>6	73.2	17.6	9.2	
Stoma depth	Flat	70.0	17.6	12.4	p = .001
	Deep	66.2	18.6	15.2	
	Very Deep	55.7	27.1	17.2	
Neck distension	No/minor	45.5	47.6	6.9	p = .00003
	Moderate/much	60.3	22.0	17.3	
Perspiration	Never/rarely	67.0	18.8	14.2	p = .0007
	Sometimes	66.5	21.1	12.5	
	Often/always	62.3	13.0	24.7	

The details provided in table 7 indicate that peristomal skin irritation is more frequent in younger patients, in earlier postoperative stages, in patients with deeper set stomas, in patients with neck distension during TE speech and in patients who experience perspiration more often. Whilst adhesive problems and skin irritation are strongly related (see Table 3 and Figure 2), the insights presented in Table 7 can help identify patients that may be more likely to develop skin irritation and provide them with support. Surgical techniques aimed at preventing deeper set stomas should be considered.

Interestingly, previous radiotherapy or chemoradiotherapy did not have an influence on reported skin irritation; in all subgroups, around two-thirds reported to never/rarely have skin irritation, 19% sometimes and 14% often/always.

2.6 Peristomal adhesives and voice rehabilitation

Sixty-five percent of the respondents (n=1844) used a VP and nearly half of them (48%) reported that the adhesive influenced their speaking, highlighting the importance of the adhesive for voice rehabilitation. TE speakers put higher demands on their adhesive due to the requirement of an airtight seal for efficient TE speech, which is reflected in the proportion that experiences problems with their adhesive: 38.5% of the VP users reported problems with the adhesives, compared to 31.2% of the non-VP users. About two-thirds of the VP users (n=1262, 68%) had tried a handsfree speaking valve, but 65% of them (n=822) had not been successful for various reasons. Problems with the adhesive accounted for nearly half (n=405, 49%) of the handsfree speaking failures. Another problem seen in patients using TE speech is 'neck-distension' during speech. Results of this research showed that 46.5% of patients with neck distension during TE speech experienced problems with their adhesive versus 30.5% of those without, indicating that neck distension negatively impacts success with peristomal adhesives.

TE speakers put higher demands on their adhesive due to the requirement of an airtight seal for efficient TE speech

3. Medical adhesives technology for TL patients

3.1 General composition of medical adhesives

Medical adhesives are generally composed of 3 layers: A carrier, also called backing, to which the adhesive ('glue') is applied, and a liner which is removed prior to application (see Figure 3). The carrier can consist of a variety of materials such as cloth, paper, or plastic. The adhesive can consist of a variety of materials such as acrylates, hydrocolloids, hydrogels, silicones, and polyurethanes. The liner is typically paper or plastic.




FIGURE 3. Layered composition of medical adhesives

The primary purpose for medical adhesives for TL patients is to place an HME in front of the tracheostoma. Secondary purposes are for example placement of a showering device or laryngectomy tube. Therefore, a fourth component is added to peristomal adhesives for TL patients: a connector for these devices, also referred to as an adaptor. The characteristics of this connector play a role in comfort and fit, as well as security of the connection with the devices it holds. Figure 4 shows the layered composition of a Provox Life Adhesive including the connector.




Figure 4. Layered composition of Provox® Life™ adhesive, including the SecureFit™ connector.

Medical adhesives are pressure sensitive, meaning that pressure applied to the surface of the adhesive will activate the adhesive by increasing the surface area contact with the skin. Over time, as the adhesive (i.e. glue) warms up, it thins and flows to fill in the gaps between the adhesive and the irregularities in the skin surface, and this increases the bond between

the adhesive and skin. The speed at which this process occurs is different for different adhesive materials. Adhesive materials such as hydrocolloids benefit from being pre-warmed, for example by rubbing them between the hands, prior to application ⁽²⁷⁾.

Medical adhesives are pressure sensitive, meaning that pressure applied to the surface of the adhesive will activate the adhesive

3.2 Material characteristics

Key factors to consider for medical adhesive materials in TL patients are:

Flexibility – more flexible and thinner materials of the medical adhesive and the adaptor can be fitted closer to the stoma and feel more comfortable in the neck. However, more rigid and thicker materials may be required to counteract neck-distension during (handsfree) TE speech.

Permeability – materials with higher permeability allow the skin to breathe and moisture to evaporate. Conversely, materials with lower or no permeability (occlusive tapes) do not allow the skin to breathe, and moisture remains trapped underneath the adhesive. In TE speakers, high permeability materials can cause undesirable air leakage through the material during speech and are therefore less suitable. This is especially the case in hands-free speech where there is no manual counterpressure keeping the adhesive materials pressed against the surface of the skin.

Absorption – hydrocolloids and hydrogels are absorbent wound dressing materials. They adhere initially to dried surfaces and as they absorb moisture adhesion weakens according to water content. These materials are considered more skin-friendly than, for example, acrylics.

Adhesion – Adhesives with lower tackiness have lower peel-strength and are therefore easier to remove from the skin. However, TL patients using TE speech often require higher tackiness for the adhesive to withstand speaking pressure.

HYDROGEL SHEET DRESSINGS

Hydrogel sheet dressings are water- or glycerin-based products, best suited for dry wounds or those with minimal to moderate exudates^(28, 29). They enable monitoring skin and/or wound without removing the dressing, as they usually are clear and translucent^(27,29,30). Furthermore, hydrogels maintain an adequate moisture balance by absorbing exudates and bacteria, and by rehydrating dry/irradiated skin⁽³¹⁾. These dressings are found to enhance epithelialization⁽³²⁻³⁴⁾, and might decrease wound infections^(35, 36). Sheet hydrogels are effective, comfortable, and soothing and create a moist wound healing environment^(27-29, 32, 37-39). The development of a hydrogel adhesive for TL patients (Provox Luna) has demonstrated benefits for overnight peristomal skin recovery and HME compliance⁽²⁰⁾. Patients using a hydrogel adhesive overnight (Provox Luna) were found to have less dermatological problems (14.3% versus 46.4% in the control group) and lower HME abandonment rates (10.7% versus 46.4% in the control group)⁽⁴⁰⁾.

HYDROCOLLOID SHEET DRESSINGS

Hydrocolloid wound dressings are occlusive or semi-occlusive dressings made of gelatin, pectin, polysaccharides, or sodium carboxymethylcellulose (CMC). The outside of the dressing is covered with a polyurethane foam or film which enables the exchange of water vapor and protects the wound against contamination from the outside^(27,41,42). The hydrocolloid dressing absorbs the wound fluid and as a result changes into a jelly-like mass, providing a moist wound healing environment and promoting autolytic debridement. Hydrocolloids are impermeable dressings with multiple benefits, including absorption of excess fluids, providing a barrier to infection, and facilitating debridement^(27,42,43).

Flexibility, permeability, absorption and adhesion are important adhesive material characteristics

MEDICAL GRADE ACRYLIC ADHESIVES

Pressure-sensitive acrylic adhesives for application to skin are made from 2-ethylhexyl acrylate, isooctyl acrylate or n-butyl acrylate copolymerized with polar functional monomers such as acrylic acid, methacrylic acid, vinyl acetate, methyl acrylate, N-vinylcaprolactam, or hydroxyethyl methacrylate. Functional comonomers increase cohesive strength, provide surface polarity, and enhance wear performance. Tack, adhesion to skin, adhesive transfer to skin, and wear performance of the adhesive are governed by the molecular weight, glass transition temperature, and the viscoelastic behavior of the adhesive. Sweating skin, a moist environment, and physical activity are the most important factors influencing the failure of an adhesive tape during wear^(27,44).

3.3 Design characteristics

Key design characteristics to consider for medical adhesives for TL patients are:

Shape of the adhesive – many tracheostomas are somewhat deep or in some cases even very deep-set into the neck. An adhesive shape that can accommodate this is beneficial and can reduce wrinkles and air tunnels which is beneficial for the effectiveness of heat and moisture exchange and TE speech.

Size of the adhesive – due to large individual variations in peristomal topography and skin condition, patient preference for size may differ. A smaller size offers a smaller contact surface area addressing possible skin concerns whereas a larger size may offer more support during speech.

Flexibility of the connector – more flexible material allows for application closer to the stoma, reducing the risk for air leaks and mucus trapping and improving comfort.

Height of the connector – a lower height reduces the profile in the neck.

Security of the connector – secure connection of the device it holds is important to prevent inadvertent cough-outs.

4. The normal skin

The skin is the largest organ in the body ⁽⁴⁵⁾ and it renews itself every 28 days ⁽⁴⁶⁾. The thickest skin is found on the soles of the feet, measuring 1.5 mm thick ⁽⁴⁷⁾.

The functions of the skin are protection (trauma, UV light, temperature, dehydration, toxins, microorganisms), maintenance of body temperature (circulatory mechanisms, sweating, insulation), sensation (pain, pressure, vibration, temperature), metabolism (produces melanin, vitamin D, excretes waste e.g. sweat) and communication (facial expression, emotions, touch) ^(47, 48).

The skin provides a semi-permeable barrier that protects the body from the external environment and maintains homeostasis in the body ^(45, 48).

The skin is the largest organ in the body

4.1 The structure of the skin

The skin consists of 3 layers (Figure 5):

- The epidermis – traps moisture ⁽⁴⁸⁾
- The dermis – provides functionality due to variety of structures in this layer ⁽⁴⁷⁾
- The subcutaneous tissue provides insulation and cushion ⁽⁴⁵⁾




FIGURE 5. The three layers of the skin: Epidermis, dermis and subcutaneous tissue.

4.1.1 The epidermis

The epidermis is the outermost visible part of the skin (Figure 6). It has a barrier function, controls water loss, and protects against UV light, bacteria and allergens ⁽⁴⁸⁾. It has a thin layer of cells varying in thickness from 0.05 mm on the eyelids to 1.5 mm on the palms and soles ⁽⁴⁷⁾.




FIGURE 6. The epidermis and its four layers.

The epidermis is the skin layer that interacts directly with peristomal adhesives and therefore it is important to understand its composition and function. It is composed of 3 types of multi-layered epithelial cells: keratinocytes, melanocytes and Langerhans cells (Figure 7). The keratinocytes account for 95% of the cells and retain moisture to maintain the skin's barrier function ⁽⁴⁵⁾, the melanocytes are responsible for skin pigmentation and protection against UV damage ⁽⁴⁵⁾, and the Langerhans cells are part of the skin's immune system ^(45, 49).

The epidermis is the skin layer that interacts directly with peristomal adhesives




FIGURE 7. The epidermis and its multi-layered epithelial cells: keratinocytes, melanocytes and Langerhans cells.

The keratinocytes layer consists of 4 types of differentiated cells (Figure 8). The innermost layer is the stratum basale (also called germinative layer) and this is where the epithelial cells (keratinocytes)

originate ⁽⁴⁸⁾. From there, the keratinocytes divide and migrate upwards through the stratum spinosum and stratum granulosum, which are living cells. As the keratinocytes migrate towards the surface they degenerate and die (and are now called corneocytes) ⁽⁴⁵⁾. The outermost layer is composed of these dead keratinocytes (corneocytes) and is called stratum corneum or sometimes the horny layer ^(49, 50). The stratum corneum is densely packed with 25 to 30 rows of flat corneocytes that give the skin structure and resilience towards bacteria and allergens ^(48, 49). The corneocytes are connected by desmosomes and they are surrounded by an extra-cellular matrix that is hydrophobic and thereby waterproofs the epidermis and prevents dehydration ⁽⁴⁸⁾. The extracellular matrix and the desmosomes allow the skin to be flexible and less susceptible to sheering forces ^(48, 51).




FIGURE 8. The four layers of differentiated keratinocytes in the outermost keratinocyte layer of the epidermis.

The outermost layer, the stratum corneum, has a 'bricks-and-mortar' like structure which helps to protect us from environmental hazards (Figure 9). The 'bricks' represent flattened, dead, keratinocytes that have matured into corneocytes and the 'mortar' represents the surrounding intercellular lipid membranes ⁽⁴⁸⁾. The bricks are tightly connected by desmosomes in overlapping layers to retain moisture while keeping allergens, pathogens and environmental toxins (e.g. UV radiation) out ⁽⁴⁸⁾. The stratum corneum is renewed about every 15 days, as dead epidermal cells are continuously shed (desquamation) ⁽⁴⁸⁾.




FIGURE 9. Graphic representation of the 'Bricks-and-Mortar' structure of the healthy epidermis

KEY TAKEAWAYS ⁽⁴⁵⁾

Epidermis must be kept intact.

Broken epidermis can allow interstitial liquid in the underlying dermis to weep to the surface. Weeping compromises application and adherence of peristomal adhesives.

If the epidermis is broken, nerve endings are exposed and will trigger pain sensation.

Pain, itching and burning sensations are signals of erosion of the epidermis and the need for prompt action.

4.1.2 The dermis

The dermis is the middle layer of skin and is about 15 to 40 times thicker than the epidermis. It consists of connective tissue, a semi-fluid substance mainly consisting of collagen and elastin fibers, and it contains nerve endings, blood vessels and skin appendages such as hair follicles, sebaceous glands and sweat glands ^(45, 53). See Figure 10.




FIGURE 10. The dermis with blood vessels, nerve endings, hair follicles, sebaceous glands and sweat glands.

The skin is a sensory organ that collects information via the nerve endings in the dermis that provide vital information about touch, pressure, pain, vibration and temperature ⁽⁴⁵⁾. The location of the sensory nerve endings in the dermis explains why superficial wounds are often more painful than wounds with more tissue loss ⁽⁵⁴⁾. See Figure 11.




FIGURE 11. Superficial location of sensory nerve endings in the dermis.

The sebaceous glands are found around the hair follicles (Figure 12) and secrete sebum, a lipid rich, oily substance which is important to maintain the moisture of the epidermis and the pH of the skin ⁽⁴⁵⁾. The secretion of sebum is controlled predominantly by sex hormones and changes with age. It is acidic and when combined with lipids it forms a hydrophobic layer that prevents skin dehydration. This helps to prevent invasion and infection by microorganisms ⁽⁴⁵⁾.




Figure 12. Sebaceous gland at hair follicle.

The sebum and sweat create an acidic environment on the skin surface in the pH range of 4 to 6 (acid), which is referred to as the acid mantle ^(45, 55) and therefore skin heals more rapidly in an acidic environment ⁽⁵⁶⁾.

KEY PERSPECTIVES ON SKIN AND MEDICAL ADHESIVES

Weeping and/or bleeding from the dermis (broken epidermis) will compromise adhesion to skin surface ⁽⁴⁵⁾.

If blood vessels are ruptured bleeding may be visible on the surface and impact adhesion ⁽⁴⁵⁾.

The dermis and epidermis are normally firmly bonded by strong collagen anchoring structures. These structures may be weakened by immune suppressive therapy. The results can be a separation of the layers when removing an adhesive, also called 'tearing'. Choice of adhesive type and use of adhesive removers are extremely important for these user groups ⁽⁵⁵⁾.

Heavy/sudden perspiration can challenge application, adhesion, and wear time. Patients

with high perspiration levels must use an adhesive with sufficient absorption capacity (i.e. hydrocolloid). When absorption capacity is reached, the adhesive will start to disintegrate and loosen. Furthermore, a more frequent changing routine may be required, e.g. after sports ⁽⁴⁵⁾.

Sebum production is important to keep the acid mantle intact and reduces with age. Choose and handle adhesive carefully and consider the age of the patient ⁽⁴⁵⁾.

Pulling hairs when removing the adhesive can cause folliculitis. Therefore, it is recommended to shave or cut hairs prior to adhesive application ⁽⁵⁷⁾.

If proper adhesive removal technique is not used, superficial layers of the skin are removed along with the adhesive product, which not only affects skin integrity but can cause pain and the risk of infection, increase wound size, and delay healing, all of which reduce patients' quality of life ⁽⁵²⁾.

4.1.3 Subcutaneous tissue

The subcutaneous tissue is the lower layer of the skin between the dermis and fascia and is composed of connective tissue, fatty tissue and contains major blood vessels and nerves (see Figure 13) ⁽⁴⁵⁾. Subcutaneous fatty tissue can often be visible in deep wounds such as pressure ulcers or traumatic wounds ⁽⁵⁸⁾. Its functions include insulation, cushion and moisture retention and its thickness depends on several factors, including the body site and age ⁽⁴⁵⁾.




FIGURE 13. Subcutaneous tissue.

5. Medical adhesive-related skin injuries (MARSI)

5.1 Categories of MARSI

Injuries related to the application or removal of adhesives are divided into the following categories: mechanical, dermatitis and other ⁽²²⁾. In Table 8 a brief description of each type of injury is provided.

TABLE 8. Categories for medical adhesive-related skin injuries. Adopted from McNichol et al. ⁽²²⁾ and Fumarola et al. ⁽⁵⁹⁾.

Mechanical
Skin stripping: Trauma that can occur after removing an adhesive incorrectly. Skin stripping removes one or more layers of the stratum corneum. Lesions are typically shallow and irregular, skin may appear shiny, and lesions may be red (erythema) or show blister formation ^(52, 60, 61) . See Figure 14 A.
Tension injury or blister: Skin injury caused by shear force as a result of stretching or movement of the skin in the presence of a solid, rigid, unbending, non-stretching adhesive, or because of inappropriate application (strapping) of the adhesive ^(60, 62, 63) .
Skin tear: An acute wound caused by trauma, cuts or scrapes, or skin friction that may become chronic. They can be partial thickness in which the epidermis separates from the dermis or full thickness in which the epidermis and dermis separate from the hypodermis ⁽⁶⁴⁾ . See Figure 14B.
Dermatitis
Irritant contact dermatitis: inflammatory skin response to chemical irritants that can cause dryness, scaling, cracking, redness, swelling, and increased skin permeability. The affected area is well-defined and corresponds with the area of exposure. The condition disturbs the skin's pH balance, which normally protects the skin from pathogens and damage by certain chemical compounds. Irritant contact dermatitis is non-allergic and commonly resolves shortly after the chemical exposure is cleared. However, the reaction may persist for weeks after exposure to the irritant. ^(65, 66) . See Figure 14C.
Allergic dermatitis: An immunological reaction in the area exposed to an allergen, for example, an adhesive component. Allergic contact dermatitis causes vesiculation (small blisters) an inflammatory response (redness, itchiness) in the area exposed and beyond, which, in turn, increases skin permeability. This type of skin irritation may be associated with other underlying allergies ^(65, 67) . See Figure 14D.
Other
Maceration: Softening and wrinkling of the skin after prolonged exposure to moisture, such as accumulation of mucus or prolonged entrapment of moisture under the adhesive. Maceration increases skin permeability and makes the skin more susceptible to infection and damage from friction and irritants. Prolonged exposure to moisture also breaks down the affected skin and loosens the adhesive. Maceration can be painful and impair wound healing. See Figure 14E.
Folliculitis: An inflammatory reaction on the skin around the hair follicles caused by a bacterial infection. Small elevations of the skin around the hair follicle that may contain pus. See Figure 14F.

The most common medical adhesive-related skin injuries in total laryngectomy patients are skin stripping, irritant contact dermatitis, allergic dermatitis and maceration ⁽²⁶⁾.



FIGURE 14 A-F. Laryngectomy stoma with skin stripping after forceful and rapid removal of the adhesive (A); example of skin tear (B); laryngectomy stoma with irritant contact dermatitis with a well-defined area of exposure to the adhesive, area is swollen and this condition is likely very painful (C); example of allergic dermatitis (D); example of maceration and desquamation on the sole of the foot (E); and example of folliculitis with small elevations around the hair follicle that contain pus (F).

5.2 Risk factors, causes and prevention

5.2.1 Intrinsic and extrinsic risk factors

Several risk factors have been identified for MARSI in the literature. They are typically divided into intrinsic and extrinsic risk factors. The factors listed in this section may not provide a complete overview and new research may uncover additional risk factors. It is therefore important to always consider the patient's general health status and any treatments they have received and medications they are taking. Knowledge of these risk factors helps identify patients that are more likely to develop MARSI.

Intrinsic risk factors:

- **Age:** Skin changes with aging include for example epidermal thinning, loss of moisture, and loss of elasticity which makes older adults more susceptible to skin injury ^(52, 68).
- **Race/ethnicity:** differences in skin structure and function among different race and ethnic groups may cause a difference in susceptibility to skin injury ⁽⁶⁹⁾.
- **(Pre-existing) Dermatological conditions:** eczema, dermatitis, chronic exudative ulcers, epidermolysis bullosa ^(22, 59).
- **Underlying medical conditions:** diabetes, infection, renal insufficiency, venous hypertension, peristomal varices ^(22, 59, 70, 71).
- **History of contact dermatitis** ⁽⁷²⁾.
- **Late post-surgical ambulation** ⁽⁷²⁾.
- **Prolonged hospital stay** ⁽⁷³⁾.
- **Low Braden scale score** ⁽⁷³⁾.

- Dry skin ⁽⁷³⁾.
- Hypoalbuminemia ⁽⁷³⁾.
- Chronic venous insufficiency ⁽⁷⁴⁾.
- Moderate to severe edema ^(59, 75).
- Hyperthermia ^(59, 75).
- Malnutrition ^(22, 59).
- Dehydration ^(22, 59).

Extrinsic risk factors:

- Drying of the skin (for example due to harsh (non pH-balanced) cleansers or low humidity) ^(22, 59).
- Prolonged exposure to moisture ^(22, 59).
- Certain medications (for example chemotherapeutic agents, long-term use of corticosteroids, anti-inflammatory agents and anticoagulants) ^(22, 59).
- Immunosuppression ⁽⁷⁶⁾.
- Radiation therapy ^(59, 77).
- Photodamage or exposure to ultraviolet light ^(22, 59).
- Repetitive adhesive tape removal ^(22, 73).
- Repeated adhesive tape application ⁽²²⁾.

Due to the nature of their underlying condition and their typically more advanced age, many TL patients may have medical conditions such as diabetes, may suffer from malnutrition or dehydration and will most likely have been exposed to radiotherapy and/or chemotherapy ⁽⁷⁸⁻⁸⁰⁾. All these factors can negatively impact their peristomal skin condition and susceptibility to MARSI.

5.2.2 Preventable causes of MARSI

Aside from the risk factors listed in the section above, a number of preventable causes of MARSI have been identified ^(22, 81, 82). These include improper choice of adhesive tape (i.e. tape with excessive adhesion for purpose, tape with insufficient stretch), improper application technique (tension ('strapping') on application, applying in wrong direction, applying to wet/moist skin, use of irritating skin barriers, which are drying to the skin, not allowing skin barriers to dry, not trimming hair, excessive use of tackifiers or glues, leaving occlusive adhesive tapes or dressings on too long, improper removal technique (i.e. quick removal, removal at high angle, insufficient support of the skin at the peel line, lack of use of adhesive removers when indicated), repeated adhesive tape/dressing changes.

5.2.3 Medical adhesive-related challenges specific for TL patients

Most of the clinical evidence around peristomal skin complications and medical adhesive-related skin injuries stems from clinical research in adjacent clinical areas such as colostomy or urostomy. Whereas these learnings are often also applicable to TL patients, there are some specific differences that deserve attention:

Speaking: Many TL patients speak with the use of a VP which requires occlusion of the HME, or the use of a handsfree speaking valve, to direct pulmonary air into the esophagus and pharynx. An airtight seal of the

evidence & experience™

adhesive is an important requirement for effective TE speech. This has the following implications for adhesive selection and application:

- Higher levels of adhesion (tackiness) are required to withstand the pressure generated when speaking and application technique should involve prevention of wrinkles.
- Material flexibility and shape of the adhesive should be considered to achieve a good fit with the peristomal skin area and minimize wrinkles and air tunnels.
- Lower moisture vapor transmission rates (permeability) of adhesive materials are required to prevent/reduce 'blow-through' during speech.
- Patients using TE speech may experience neck distension during speech, especially when using a handsfree speaking valve. An adhesive with a more firm/rigid adaptor may be helpful in supporting the neck tissues and limiting the extent of the distension.

Mucus: TL patients cough up mucus (Airway Surface Liquid – ASL) via their tracheostoma. This is ideally cleared away immediately, but in some instances, mucus may accumulate and/or seep underneath the adhesive see Figure 15. The pH of ASL averages about 6.6 and can be more alkaline in patients experiencing pulmonary disease ⁽⁸³⁾. The prolonged presence of ASL may affect the acid mantle of the peristomal skin which may cause skin maceration and skin irritation and it negatively affects the airtight seal of the adhesive



FIGURE 15. Mucus accumulation underneath the adhesive causing the adhesive to come off and possibly causing maceration

Movement: the tracheostoma is located at the base of the neck and therefore the adhesive baseplate can be affected by movements of the head and neck. Patients may experience discomfort from the adhesive during neck movement.

Appearance: The tracheostoma at the base of the neck is in a visible location and hence the appearance of devices used in the area is of importance.

6. Provox® Life™ Adhesives

Provox Life adhesives were designed and developed taking into consideration the learnings from the ethnographic and market research conducted by Atos Medical, scientific literature, and utilizing the most recent and advanced adhesive material technologies. Adhesives were developed using 3 materials: hydrocolloid, hydrogel and acrylic and in a novel clover-shape. Furthermore, the Provox Life Adhesives are equipped with patented SecureFit connector technology, and one adhesive is designed in a concave shape for deeper set tracheostomas, with a firm base and stabilizing bars to provide support to the neck during (handsfree) TE speech.

6.1 Material characteristics and technical details

In Table 9 the material characteristics of the Provox Life Adhesives are detailed, alongside with some characteristics of the type of material used in general. Please refer to section 3.2 for details on material characteristics of hydrocolloids, hydrogels and medical acrylics in general.

TABLE 9. Technical details Provox Life Adhesives (data on file), including general characteristics of the types of adhesive materials used (see Section 3.2) and relevant differences with the legacy Provox Adhesives (data on file).

Provox Life Adhesive	Technical details
Sensitive Adhesive Skin-friendly adhesive for everyday use  	Carrier: Non-perforated ethyl methyl acetate. Adhesive: Hydrocolloid. Synthetic rubber base with a proprietary mixture of: Tackifiers - to create the (initial) adhesion to the skin surface Antioxidants - to prevent ageing of the material due to UV light and temperature Absorbents - to absorb moisture (this is the specific component called hydrocolloid) Liner: Paper with siliconized surface. One piece. Connector: SecureFit. Diameter 23 mm. Application: The adhesive material (hydrocolloid) is pressure-sensitive and reaches maximum adhesion after about 15 minutes. Hydrocolloid material responds well to pre-warming prior to application which helps increase initial adhesion to the skin surface. Removal: Once the hydrocolloid material absorbs moisture, it becomes a gel and loosens from the skin surface, making it easier to remove. Use Provox Adhesive Remover wipe underneath adhesive – starting from finger lift tab. Or use a lint-free cloth with water (distilled, sterile or drinking quality), apply at the edges underneath the adhesive and allow time for the adhesive to absorb the water and loosen. Technical data versus legacy device Provox OptiDerm: The material is 28% thinner, which increases flexibility and improves comfort and fit. Material has a higher level of integration which makes it more robust towards break down from moisture. Sizes: Comes in 3 sizes- Round, Oval, Plus.
Night Adhesive Soothing adhesive for skin rest and recovery  	Carrier: Polyurethane. Non-perforated. Adhesive: Hydrogel with integrated polypropylene non-woven reinforcement layer to maintain integrity of the adhesive. Soft and high tack, designed especially for skin contact with and for use on breached or damaged skin. The hydro-active adhesive has the following characteristics: Highly breathable – hydrogels have significantly more breathability than other materials such as silicones, acrylics, or hydrocolloids owing to a high moisture transport rate (MTR), which is a measure of breathability. Moisturizing - the hydrogel is formulated with ~30% pharmaceutical grade, vegan glycerin which is designed to maintain skin moisture levels as well as the integrity of the adhesive, making it less prone to breaking or tearing. Water soluble - hydrogels are water soluble and can be rinsed or soaked off the skin. This makes the adhesive easy to remove, but also less suitable for prolonged exposure to external fluids (e.g., showering). Wound cooling effect – the high-water content of hydrogels provides a cooling effect. Absorption of bacteria - hydrogels are capable of absorbing bacteria within their polymer structure which may help prevent the accumulation of excess bacteria within the top layers of the skin. Transparent – hydrogel materials are transparent which allows for observation of potential wound or skin defects underneath the adhesive and is discreet. Thin and flexible – the adhesive is thin and flexible which adds to patient comfort. Liner: Polyester with siliconized surface. Three-piece. Connector: SecureFit. Diameter 23 mm. Packaging: High moisture barrier packaging to preserve water content. Application: Pressure sensitive material. Does not benefit from pre-warming prior to application. Must be applied to clean and dry skin, free of skin barriers, creams, oily residues, etc. Removal: Use lint-free cloth with saline or water (distilled, sterile or drinking quality) to remove. Technical data versus legacy device Provox Luna Adhesive:The materials used are the same. The diameter of the connector is 23 mm instead of 19 mm and is therefore compatible with all Provox Life HMEs. Sizes: Available in one size only

Standard Adhesive Flexible adhesive with strong adhesion  	Carrier: Polyethylene. Perforated. Adhesive: Medical Grade Acrylic Liner: Paper with siliconized surface. One piece. Connector: SecureFit. Diameter 23 mm. Application: The pressure sensitive medical grade acrylic material reaches maximum adhesion after about one hour and does not benefit from pre-warming prior to application. Removal: Reduced MTR may lead to higher levels of moisture trapping underneath the adhesive which helps to reduce stickiness for adhesive removal. Use Provox Adhesive Remover wipe on top of perforated carrier to saturate material and loosen adhesive prior to removal. Then use adhesive remover wipe from the side to gently remove adhesive - starting from finger lift tab. Technical data versus legacy device Provox Flexiderm adhesive: The materials of this adhesive are the same as for the legacy Provox FlexiDerm device. However, the coat weight (amount of adhesive (i.e., glue)) is 67% higher which increases stickiness, lowers permeability, and reduces MTR. Lower permeability and MTR reduces blow-through (i.e., air blowing through the adhesive material in spots where it is not connecting with the peristomal skin) during TE speech. Lab-tests demonstrated a >400% reduction in blow-through. Sizes: Comes in 3 sizes- Round, Oval, Plus.
Stability Adhesive Firm, concave adhesive providing support and strong adhesion  	Carrier: Polyethylene with siliconized surface. Perforated. Adhesive: See above. Same as Provox Life Standard Adhesive. Liner: Polypropylene. Three-piece. Connector: SecureFit. Diameter 23 mm. Firm, concave base with stabilizing bars for deeper set stomas and neck support during (handsfree) TE speech. Application and Removal: See above. Same as Provox Life Standard Adhesive. Technical data versus legacy device Provox StabiliBase: The design of the firm, concave base and connector has remained the same. For materials, see information in section above for Provox Standard). Sizes: Available in one size only

6.2 Shape

The flat Provox Life Adhesives (i.e., Sensitive Adhesive and Standard Adhesive) come in 3 sizes (Round, Oval and Plus) and Night in 1 size (Oval). The variation in sizes helps to address differences in stoma configurations and patients’ needs and preferences. All are clover-shaped, in which each leaf of the clover moves independently to improve fit to the stoma and reduce air tunnels and wrinkles in the adhesive. See Figure 16.




FIGURE 16. Simulation of improved fit of clover-shaped design with independent movement of the leaves.

6.3 Connector

The Provox Life SecureFit connector is soft and low-profile and the patented SecureFit technology provides a secure connection for the devices it holds in front of the tracheostoma. The soft and flexible material of the connector is designed to be comfortable and allows application of the adhesive close to the stoma which helps to obtain a tight seal around the stoma for a more secure fit. Whilst the inner diameter of the connector is wider than the legacy Provox Adhesives, the outer diameter of the connector is smaller which allows application of the adhesive closer to the stoma. See Figure 17.




FIGURE 17. Dimensions of the Provox connector and the Provox Life SecureFit connector.

7. Clinical evidence

Clinical data on Provox Life is available from one large user experience program conducted by Atos Medical (data on file) and two Post-Market Clinical Follow Up studies sponsored by Atos Medical ^(4, 5). The sections below provide details about the clinical data regarding the Provox Life Adhesives.

7.1 Post-market User Experience Program

During 2020 and 2021 Atos Medical conducted a large scale, structure Post Market Experience Program for Provox Life that included 437 patients in 10 countries (Denmark, Finland, Sweden, the Netherlands, Germany, Switzerland, UK, US, New Zealand, Italy). The average age was 65.4 years (range 36-89), average postop follow up was 5 years (range 0 – 27 years) and the male-female ratio was 8.1 – 1.9. The timing and length of follow up varied slightly between countries and due to weather-related supply and follow up challenges the follow-up response rate in the US was low. Out of the 437 patients approached, 364 patients completed the baseline questionnaire, 187 the 2-week follow up, 180 the 6-week follow up and 155 the 12-week follow up questionnaire.

Results showed a significant impact of Provox Life on pulmonary health, as demonstrated by significant reductions in Coughing and Sputum symptoms and impact as per the validated Cough And Sputum Assessment Questionnaire (CASA-Q) ⁽¹⁴⁾.

With regards to the HME attachments, the percentage of patients using adhesives increased from 76.4% for the usual routine, to 81.1% at 2 weeks, 82.5% at 6 weeks and 86% at 12 weeks. The device life of the Provox Life Adhesives did not differ significantly from the device life during usual routine, with averages of 29.1 hours for usual care (legacy Provox Adhesive), 29.1 hours at 2 weeks, 29.1 hours at 6 weeks and 27.4 hours at 12 weeks. Notably, the median device life remained 24 hours throughout, indicating that most patients adhere to a daily routine. The reasons to change the adhesive fluctuated over time and showed an overall increase in routine changes. Of note, as the percentage of routine changes increases, those of replacement for leakage and blow-through decrease (see Figure 18). Patient preference and education may dictate whether the adhesive is changed when it fails versus routinely on a set schedule.

In addition to the Provox Life Adhesives, patients also received skin care products to support peristomal skin health. Patients were educated to use Provox Cleaning Towel to clean the skin prior to application and Provox Skin Barrier wipes to apply a protective layer over the skin prior to applying the adhesive.




FIGURE 18. Shifts in adhesive change patterns from usual care and over time with the Provox Life Adhesives during the Post-market Experience Program conducted by Atos Medical (data on file).

Provox Adhesive Remover wipes were recommended for careful and gentle removal of the acrylic and hydrocolloid adhesives and water was recommended for the gentle and careful removal of hydrogel adhesives. The graph in Figure 19 indicates an increase in the use of peristomal skin care products.




Figure 19. Usage of peristomal skin care products over time when patient switched to Provox Life during the Post-market Experience Program conducted by Atos Medical (data on file).

Using Provox Life Adhesives, the frequency of skin irritation reduced significantly with 49.7% reporting sometimes/often/always skin irritation during Usual Routine versus 32.4% at 12 weeks for Provox Life (Friedman paired test, with Bonferroni post hoc correction for multiple analyses, $P<.001$). See Figure 20.




FIGURE 20. Reduced skin irritation following the change from Usual Routine to Provox Life Adhesives during Post-market Experience Program conducted by Atos Medical (data on file).

After 12 weeks, 88% of the patients preferred the Provox Life Adhesives over their usual adhesive.

Post-market surveillance data showed a reduction in the frequency of skin irritation with Provox Life Adhesives

7.2 Post-market clinical follow-up studies

The first clinical study with Provox Life was conducted in Italy by Longobardi et al ⁽⁴⁾. The study was a cross-over clinical trial including 40 laryngectomized patients with random allocation to the sequence. Group 1 started with 6 weeks of Provox Life followed by 6 weeks of Usual Care and group 2 started with 6 weeks of Usual Care followed by Provox Life. Results with regards to the adhesives showed that the average device life was not significantly different (20.3 hours for Usual Care, 19.3 for Provox Life). Frequency of skin irritation reduced significantly with Provox Life, with 35% of participants reporting never/rarely experiencing skin irritation during Usual Care and 58% during Provox Life. This data was confirmed by a diary, in which subjects reported skin irritation on average 4.3 days during the last 2 weeks of the UC period and 2.6 days during the last 2 weeks of using Provox Life devices (31% vs 19% of the days, $P = .013$).

Number of days with skin irritation reduced significantly with Provox Life Adhesives

A second clinical study was conducted in Australia by Ward et al ⁽⁵⁾. This study was a prospective, multi-center two-phase study that investigated the impact of implementing an optimal day-night routine on pulmonary health. Forty-two patients completed the study. During phase 1, patients changed from their usual HME and attachment to a similar Provox Life HME and attachment, during Phase 2 the optimal day-night routine was implemented, and a wider variety of devices was used. There was no significant difference in device life of the adhesives between usual routine (38.3 hours), phase 1 (37.2 hours) and phase 2 (36.1 hours). The frequency of skin irritation did not differ significantly between the usual routine and the two study phases and remained at a low level throughout the study. The hydrogel adhesive was used during the night to alleviate skin irritation by 36% of the patients, with 24% using it every night and 12% every other night.

About one-third of TL patients used a hydrogel adhesive overnight to alleviate skin irritation

Longobardi et al., 2022



Ward et al., 2023



8. Prevention of MARSI in TL patients

As indicated in section 5.2.2 several preventable causes of MARSI have been identified (i.e. improper choice of adhesive, improper application technique, leaving non-permeable adhesives on too long, improper removal technique, repeated adhesive changes, touching the adhesive layer prior to application). In this section, each of these causes is addressed together with possible solutions to consider in general and for TL patients in particular. In addition, some causes for MARSI that are specific to TL patients are addressed here as well: Speaking too soon after application, mucus trapping underneath adhesive, excessive speaking pressure, not supporting the adhesive when removing the HME.

Table 10. Preventable causes of MARSI in TL patients and possible solutions.

Preventable causes to medical adhesive-related skin injuries		Recommendations and solutions to consider
Improper choice of adhesive	Excessive adhesion (tackiness) for the purpose/situation. Wrong choice of adhesive	Always select the most skin-friendly option based on the patients' skin and communication needs. Provox Life Sensitive Adhesive and Provox Life Night Adhesive are the most skin-friendly options. Provox Life Sensitive Adhesive is a skin friendly hydrocolloid that is suitable for immediate postoperative use and can be used on an ongoing basis if the tackiness is sufficient for the patients' communication needs. Provox Life Night Adhesive is a hydrogel adhesive that can be used overnight on dry skin. Hydrogel contains water and provides a moist environment for the skin to recover. This hydrogel adhesive can also be used on an ongoing basis during the day and night if the tackiness is sufficient for the patients' communication needs. It is best suited for dry skin. See section 3.2 for more details on the materials characteristics of hydrocolloids and hydrogels.
Improper application technique	Applying with tension Applying in wrong direction Applying to wet/moist skin Applying to greasy or oily skin Allowing for gaps and wrinkles Using irritating skin barriers (these are drying to the skin) Not allowing skin barriers to dry (common cause of irritant contact dermatitis) Not trimming hair Excessive use of glues or tackifiers	Apply the adhesive with the neck in neutral position, do not pull or stretch the adhesive when applying. Start from the center, around the stoma, then outwards starting with the bottom then top and sides. Ensure the skin is clean and dry, free of mucus and wound exudate. Cough and clear secretions prior to application. If the skin is wet due to wound exudate, use a moisture absorbing hydrocolloid adhesive such as Provox Life Sensitive Adhesive and ensure the skin is dry prior to adhesive application. Ensure the skin is clean and dry, free from creams, ointments, and adhesive removers. Avoid gaps and wrinkles that can allow moisture and mucus to enter between the adhesive and the skin, which can create an air leak during TE speech. Use no-sting skin barriers such as Provox Skin Barrier that are suited for peristomal skin and for use around a neck stoma. Skin barriers with strong smells of vapors and skin barriers sprays should be avoided as they can be inhaled. Let skin barriers dry fully prior to adhesive application. Trim/shave all hair in and around the area of application. Use tackifiers and silicone glues only when necessary and do so sparingly.
Leaving non-permeable adhesives on too long		Establish a change routine that optimizes for a good and secure seal for the duration of a full day, allowing for scheduled adhesive changes at a convenient time. Patients may define success by the length of time the adhesive stays in place. This can cause skin irritation due to the adhesive staying on too long, disruption due to unscheduled adhesive changes, poor application due to unanticipated adhesive changes under poor and unhygienic conditions, and insecurity due to the adhesive losing the seal at unpredictable times, disrupting the ability to communicate effectively. Scheduled routine changes allow for properly prepared adhesive changes under good and hygienic conditions.
Improper removal technique	Quick removal	Remove the adhesive slowly and gently, starting from the finger-lift tab and from the bottom of the stoma upwards

Table 10. Preventable causes of MARSI in TL patients and possible solutions.

Preventable causes to medical adhesive-related skin injuries	Recommendations and solutions to consider
Removal at high angle (i.e. perpendicular at 90°) Lack of support to skin at peel line Not using adhesive remover when indicated Not using water to remove the adhesive when indicated	Peel the adhesive off at a low angle (i.e. close to 180°, not straight up (see Figure 21). This reduces sheer forces on the skin. Support the skin at the peel line (see Figure 21) Use no-sting adhesive remover wipes free of strong smells and vapors, suitable for peristomal skin and tracheostomas, such as Provox Adhesive Remover. Application method differs based on the adhesive materials. On permeable adhesives such as the acrylic Provox Life Standard Adhesive, first apply the adhesive remover on top of the adhesive. The adhesive remover will permeate the adhesive carrier and help dissolve the glue. Then apply the adhesive remover underneath the adhesive whilst gently peeling the adhesive off the skin at a low angle. In this case, the finger with the adhesive remover wipe provides support at the peel line at the same time. Adhesive remover may be used to remove hydrocolloid adhesives such as Provox Life Sensitive Adhesive. This adhesive is non-permeable, so the adhesive remover will only be applied during removal, at the edge of the adhesive. Hydrocolloids may also be removed with water. The water must be applied at the edge of the adhesive and sufficient time must be allowed for the hydrocolloid to absorb the water. Hydrogels are removed with water. They take up water rapidly and the hydrogel will easily release from the skin when water is applied at the edge of the adhesive.
Repeated adhesive changes	Avoid repeated adhesive changes, for example multiple times per day. Identify the cause of the repeated changes and work towards a sustainable routine.
Speaking too soon after application	Medical adhesives are pressure sensitive; apply firm pressure to the surface to activate the adhesive by increasing the surface area contact. Over time, the adhesive (glue) will warm and flow to fill in gaps between the adhesive and irregularities in the skin surface which increases the strength of the bond. The length of time this takes differs for the different materials. Pre-warm hydrocolloids prior to application to help speed this process.
Mucus trapping underneath adhesive	Do not allow secretions (mucus) from coughing to settle at the bottom of the stoma and seep underneath the adhesive. After a productive cough, any secretions must be removed from the area.
Excessive speaking pressure	Excessive and sudden speaking pressure puts strain on the adhesive and can cause premature loosening of the adhesive. Work on regulating and lowering speaking pressure and provide medical treatment for hypertonicity if indicated.
Not supporting the adhesive when removing the HME	Removing the HME from the connector may pull the adhesive away from the skin. This may cause skin irritation and premature loosening of the adhesive leading to repeated adhesive changes. Always support the adhesive with two fingers (one left and one right) when removing the HME.
Touching the adhesive layer prior to application	Minimize touching the sticky adhesive layer. Finger contact may decrease the adhesion capacity of the adhesive. This may cause blow-through and premature loosening of the adhesive leading to repeated adhesive changes. For adhesives with stepwise peel-off liners, using the stepwise method will help to avoid this finger contact.

As is evident from the table above, good skin care, proper adhesive selection and appropriate application and removal routine and technique are all important components of successful adhesive use and help to avoid MARSI.




FIGURE 21. A proper adhesive removal technique includes slow and gentle removal with a low-angle peel direction and finger support at the peel line, supported with adhesive remover wipes when indicated.

9. Conclusions

Healthy peristomal skin is a prerequisite for optimal medical adhesive use. It creates the foundation for adherent HME use, which leads to improved pulmonary health and facilitates tracheoesophageal voicing (see Figure 22).




FIGURE 22. Promotion of good peristomal skin health and prevention and timely management of skin problems supports adherence to HME use and an airtight seal benefiting voice and pulmonary rehabilitation.

Minimization of Peristomal Skin Complications and Medical Adhesive-Related Injuries helps keep the skin intact.

This Whitepaper provides knowledge and tools that can help HCP promote good skin health and educate their patients on adhesive selection and proper application and removal techniques.

10. Q & A with the expert



Prof Dr Thomas Rustemeyer

Biography

Thomas Rustemeyer is professor in dermatology and allergology and occupational dermatology at the Department of Dermatology-Allergology at the Amsterdam University Medical Centers, the Netherlands. He is supervisor of immunological research of the University, and he is active in supervising patient's care and research addressing fundamental immunological and translational research topics. He is (co)author of more than 300 publications and editor of several textbooks. He is member and/or chairperson of various national and international societies active in dermatology, allergology, immunology and occupational health, including governmental and regulatory affairs.

Question:

Are some patients more at risk for skin breakdown? And if so, how can we identify them?

Answer:

This is a very good point – to prevent problems. Not everyone has the same skin condition. Skin tends to get thinner as we age. As people get older their skin becomes more fragile. And the same holds true for people who had radiotherapy at the site and people who use topical steroids that cause thinning of the skin. In addition to age, diabetes can also slow down skin regeneration and weakens the skin. Thus, people who suffer from long-term diabetes also have fragile skin and an increased risk for wound infection. In these cases, you should pay extra attention not to harm the patient by removing the adhesive with too much force or by leaving it on for too long.

Question:

How do you identify whether something is an irritation or an infection of the skin?

Answer:

This is not black and white, because infection follows irritation, and sometimes irritation can follow infection, but this is not so frequent. So, to prevent infection, you need to prevent irritation and therefore these

evidence & experience™

baseplates are so important, even using them as soon as possible after surgery.

Irritation from the adhesive should be sharply demarcated, so where the adhesive is, the irritation should end. An infection is causing a reaction of the immune system, so it tends to be broader and becomes bigger within 1 or 2 days, it's growing as it is an infection.

Question:

Do you have any tips on removing the adhesive?

Answer:

It must be done gently and carefully, so take your time. Do not remove it abruptly and quickly as this will pull off the horny cells which are the skin barrier. In addition, the sheering forces can lead to irritation of the skin. Inform the patient to take their time to remove the adhesive. If possible, in front of a mirror. Gently peel it off millimeter by millimeter. Leave the skin open for a couple of minutes so it can rest and evaporate. Maybe even a little longer, to help the skin settle down. These are the most important things.

I would recommend not to use any creams or body lotions, because they have preservatives in them, and they are often scented. Then you have the fragrances and preservatives on the skin and when you apply the adhesive over that you have a higher penetration of the allergens which in the end may lead to allergic contact dermatitis. One final thing, you could be allergic to the adhesive, but the risk is very little. If a patient tells you and starts to complain that the skin is getting worse and worse, think about the contact allergy. You can simply test it by putting a piece of adhesive, for example on the arm or the back, a patch test, and leave it there for a few days. And if the patient also gets itchy, red and eczematous skin there, where the skin was normal, also gets itchy, red and eczematous skin, then you are most likely dealing with allergic dermatitis. And this irritation can be very well treated by using corticosteroids (from a nose spray) underneath the baseplate. Use nose spray, not cream or ointment, because they are greasy and will prevent the adhesive from sticking.

Acknowledgment

We would like to thank Prof Dr Rustemeyer for reviewing the content of this Whitepaper and for his valuable comments.

Any recommendations in this educational material are a general guide for best practice, to be implemented by qualified healthcare professionals subject to clinical judgement and availability of healthcare resources.

The information presented should not be considered medical advice for specific conditions. A patient's individual circumstances and preferences should always be considered and clinical practice should be in accordance with the principles of protection, participation and partnership.

Bibliography

- Bien S, Okla S, van As-Brooks CJ, Ackerstaff AH. The effect of a Heat and Moisture Exchanger (Provox HME) on pulmonary protection after total laryngectomy: a randomized controlled study. *Eur Arch Otorhinolaryngol*. 2010;267(3):429-35.
- Parrilla C, Minni A, Bogaardt H, Macri GF, Battista M, Roukos R, et al. Pulmonary Rehabilitation After Total Laryngectomy: A Multicenter Time-Series Clinical Trial Evaluating the Provox XtraHME in HME-Naive Patients. *Ann Otol Rhinol Laryngol*. 2015;124(9):706-13.
- Beck ACC, Retel VP, Bunting G, Sethi R, Deschler DG, van den Brekel MWM, et al. Cost-effectiveness analysis of using the heat and moisture exchangers compared with alternative stoma covers in laryngectomy rehabilitation: US perspective. *Head Neck*. 2020;42(12):3720-34.
- Longobardi Y, Galli J, Di Cesare T, D'Alatri L, Settimi S, Mele D, et al. Optimizing Pulmonary Outcomes After Total Laryngectomy: Crossover Study on New Heat and Moisture Exchangers. *Otolaryngol Head Neck Surg*. 2022;167(6):929-40.
- Ward EC, Hancock K, Boxall J, Burns CL, Spurgin AL, Lehn B, et al. Post-laryngectomy pulmonary and related symptom changes following adoption of an optimal day-and-night heat and moisture exchanger (HME) regimen. *Head Neck*. 2023;45(4):939-51.
- Kuo M, Ho CM, Wei WI, Lam KH. Tracheostomal stenosis after total laryngectomy: an analysis of predisposing clinical factors. *Laryngoscope*. 1994;104(1 Pt 1):59-63.
- Verschuur HP, Gregor RT, Hilgers FJ, Balm AJ. The tracheostoma in relation to prosthetic voice rehabilitation. *Laryngoscope*. 1996;106(1 Pt 1):111-5.
- Hennessey M, Bann DV, Patel VA, Saadi R, Kreml GA, Deschler DG, et al. Commentary on the management of total laryngectomy patients during the COVID-19 pandemic. *Head Neck*. 2020;42(6):1137-43.
- Iszlai Z, Fodor B, Szabo R, Szekanecz Z, Karosi T. [Thermal imaging and velocity measurements of the exhaled airflow in total laryngectomized patients during COVID-19 pandemic.]. *Orv Hetil*. 2023;164(34):1327-36.
- Herranz J, Espino MA, Morado CO. Pulmonary rehabilitation after total laryngectomy: a randomized cross-over clinical trial comparing two different heat and moisture exchangers (HMEs). *Eur Arch Otorhinolaryngol*. 2013;270(9):2479-84.
- Dassonville O, Mérol JC, Bozec A, Swierkosz F, Santini J, Chais A, Marcy PY, Giacchero P, Chamorey E, Poissonnet G. Randomised, multi-centre study of the usefulness of the heat and moisture exchanger (Provox HME®) in laryngectomised patients. *Eur Arch Otorhinolaryngol*. 2011 Nov;268(11):1647-54..
- Hilgers FJ, Aaronson NK, Ackerstaff AH, Schouwenburg PF, van Zandwijk N. The influence of a heat and moisture exchanger (HME) on the respiratory symptoms after total laryngectomy. *Clin Otolaryngol Allied Sci*. 1991;16(2):152-6.
- Ackerstaff AH, Fuller D, Irvin M, Maccracken E, Gaziano J, Stachowiak L. Multicenter study assessing effects of heat and moisture exchanger use on respiratory symptoms and voice quality in laryngectomized individuals. *Otolaryngol Head Neck Surg*. 2003;129(6):705-12.
- Atos Medical. Whitepaper: Provox Life - Closing in on the humidification deficit in total laryngectomy patients. 2021.
- Pedemonte-Sarrias G, Villatoro-Solagaistoa JC, Ale-Inostroza P, Lopez-Vilas M, Leon-Vintro X, Quer-Agusti M. Chronic adherence to heat and moisture exchanger use in laryngectomized patients. *Acta Otorrinolaringol Esp*. 2013;64(4):247-52.
- Bassiouni M, Wakonig KM, Dommerich S. [Surgical refinements towards a functional voice prosthesis and cannula-free tracheostoma in total laryngectomy]. *Laryngorhinootologie*. 2021;100(3):174-5.
- Hilgers FJ, Ackerstaff AH, Van As CJ, Balm AJ, Van den Brekel MW, Tan IB. Development and clinical assessment of a heat and moisture exchanger with a multi-magnet automatic tracheostoma valve (Provox FreeHands HME) for vocal and pulmonary rehabilitation after total laryngectomy. *Acta Otolaryngol*. 2003;123(1):91-9.
- Lansaat L, de Kleijn BJ, Hilgers FJ, van der Laan BF, van den Brekel MW. A prospective multicenter clinical feasibility study of a new automatic speaking valve for postlaryngectomy voice rehabilitation. *Eur Arch Otorhinolaryngol*. 2016;274(2):1005-13.
- Op de Coul BM, Ackerstaff AH, van As-Brooks CJ, van den Hoogen FJ, Meeuwis CA, Manni JJ, et al. Compliance, quality of life and quantitative voice quality aspects of hands-free speech. *Acta Otolaryngol*. 2005;125(6):629-37.
- Ratnayake CBB, Fles R, Tan IB, Baijens LWJ, Pilz W, Meeuwis CA, et al. Multicenter randomized crossover trial evaluating the provox luna in laryngectomized subjects. *Laryngoscope*. 2019;129(10):2354-60.
- Brook I, Bogaardt H, van As-Brooks C. Long-term use of heat and moisture exchangers among laryngectomees: medical, social, and psychological patterns. *Ann Otol Rhinol Laryngol*. 2013;122(6):358-63.
- McNichol L, Lund C, Rosen T, Gray M. Medical adhesives and patient safety: state of the science: consensus statements for the assessment, prevention, and treatment of adhesive-related skin injuries. *J Wound Ostomy Continence Nurs*. 2013;40(4):365-80; quiz E1-2.
- Leemans M, van Sluis KE, van Son R, van den Brekel MWM. Interaction of functional and participation issues on quality of life after total laryngectomy. *Laryngoscope Investig Otolaryngol*. 2020;5(3):453-60.
- Nybaek H, Jemec GB. Skin problems in stoma patients. *J Eur Acad Dermatol Venereol*. 2010;24(3):249-57.
- Hoefloek J, Guy D, Allen S, St-Cyr D. A prospective multicenter evaluation of a moldable stoma skin barrier. *Ostomy Wound Manage*. 2009;55(5):62-9.
- Leemans M. General discussion and future perspective. Sticking to it: HME humidification, optimization and fixation for laryngectomized patients: University of Amsterdam; 2023. p. 186-7.
- Rippon M, White R, Davies P. Skin adhesives and their role in wound dressings. *Wounds UK*. 2007;3(4):76-86.
- Sarabahi S. Recent advances in topical wound care. *Indian J Plast Surg*. 2012;45(2):379-87.
- Hess CT. Addressing wounds. *Nurs Manage*. 2000;31(9):73-5.
- Hess CT. When to use hydrogel dressings. *Adv Skin Wound Care*. 2000;13(1):42.
- Gollins S, Gaffney C, Slade S, Swindell R. RCT on gentian violet versus a hydrogel dressing for radiotherapy-induced moist skin desquamation. *J Wound Care*. 2008;17(6):268-5.
- Sood A, Granick MS, Tomaselli NL. Wound Dressings and Comparative Effectiveness Data. *Adv Wound Care (New Rochelle)*. 2014;3(8):511-29.
- Beam JW. Occlusive dressings and the healing of standardized abrasions. *J Athl Train*. 2008;43(6):600-7.
- Pars H, Cavusoglu H. Effects of 3 Different Methods of Care on the Peristomal Skin Integrity of Children with Percutaneous Endoscopic Gastrostomy Tubes: A Prospective Randomized Controlled Trial. *Adv Skin Wound Care*. 2018;31(4):172-81.
- Blumenstein I, Borger D, Loitsch S, Bott C, Tessmer A, Hartmann F, et al. A glycerin hydrogel-based wound dressing prevents peristomal infections after percutaneous endoscopic gastrostomy (PEG): a prospective, randomized study. *Nutr Clin Pract*. 2012;27(3):422-5.
- Aschl G, Kirchgatterer A, Fleischer M, Hinterreiter M, Hubner D, Kranewitter W, et al. [The frequency of wound infections after PEG-placement and utilization of glycolgel wound dressing: a randomized controlled trial]. *Wien Klin Wochenschr*. 2008;120(7-8):224-7.
- Hampton S. A focus on ActiFormCool in the reduction of pain in wounds. *Br J Community Nurs*. 2007;12(9):S37-S42.
- Williams C. Hydrogel dressings: their role in dry and sloughy wounds. *Community Nurse*. 1998;4(9):42-4.
- Tam KW, Lai JH, Chen HC, Hou WH, Ko WS, Chen SL, et al. A systematic review and meta-analysis of randomized controlled trials comparing interventions for peristomal skin care. *Ostomy Wound Manage*. 2014;60(10):26-33.

40. Mayo-Yanez M, Balboa-Barreiro V, Lechien JR, Calvo-Henriquez C, Chiesa-Estomba C, Herranz Gonzalez-Botas J, et al. Improving heat and moisture exchanger therapy with a hydrogel base adhesive in laryngectomized patients: an open randomized crossover trial. *J Laryngol Otol*. 2022;1-20.
41. Rajendran S. *Advanced Textiles for Wound Care*. Elsevier Science; 2018.
42. Sood A, Granick MS, Tomaselli NL. Wound Dressings and Comparative Effectiveness Data. *Adv Wound Care (New Rochelle)*. 2014;3(8):511-29.
43. Thomas S. Hydrocolloid dressings in the management of acute wounds: a review of the literature. *International Wound Journal*. 2008;5(5):602-13.
44. Kenney JF, Haddock TH, Sun RL, Parreira HC. Medical-grade acrylic adhesives for skin contact. *Journal of Applied Polymer Science*. 1992;45(2):355-61.
45. Shimizu H. Structure and function of the skin. *Shimizu's Dermatology*. 2017. p. 1-42.
46. Shimizu H. Disorders of keratinization. *Shimizu's Dermatology*. 2017. p. 285-324.
47. Burian EA. Huden. In: Susan Bermark BØM, editor. *Sår og sårbehandling : En grundbog i sygeplejen*. 1: FADL's Forlag; 2017. p. 13-28.
48. Langøen AaB, J. Maintaining Skin Integrity. In: Flanagan M, editor. *Wound Healing and Skin Integrity: Principles and Practice*. West Sussex: John Wiley & Sons, Ltd; 2013. p. 18-32.
49. Ersser SJ, Getliffe K, Voegeli D, Regan S. A critical review of the inter-relationship between skin vulnerability and urinary incontinence and related nursing intervention. *Int J Nurs Stud*. 2005;42(7):823-35.
50. Harding CR, Watkinson A, Rawlings AV, Scott IR. Dry skin, moisturization and corneodesmolysis. *Int J Cosmet Sci*. 2000;22(1):21-52.
51. Janniger CK, Schwartz RA, Szepietowski JC, Reich A. Intertrigo and common secondary skin infections. *Am Fam Physician*. 2005;72(5):833-8.
52. Cutting KF. Impact of adhesive surgical tape and wound dressings on the skin, with reference to skin stripping. *J Wound Care*. 2008;17(4):157-8, 60-2.
53. Baroni A, Buommino E, De Gregorio V, Ruocco E, Ruocco V, Wolf R. Structure and function of the epidermis related to barrier properties. *Clin Dermatol*. 2012;30(3):257-62.
54. Saffle JR, Banister, M., Cahalane, M., Lee, J. O., Palmieri, T. L. Burns. In: Lawrence PF, editor. *Essentials of General Surgery*. 5th ed. Philadelphia: Lippincott Williams & Wilkins; 2013. p. 188.
55. Lyon CC, Smith, Amanda J. *Abdominal Stomas and Their Skin Disorders: An Atlas of Diagnosis and Management*. 1st ed: CRC Press; 2001.
56. Aschl G, Kirchgatterer A, Fleischer M, Hinterreiter M, Hubner D, Kranewitter W, et al. [The frequency of wound infections after PEG-placement and utilization of glycolgel wound dressing: a randomized controlled trial]. *Wien Klin Wochenschr*. 2008;120(7-8):224-7.
57. Cowley K, Vanoosthuyze K. The biomechanics of blade shaving. *Int J Cosmet Sci*. 2016;38 Suppl 1:17-23.
58. Beeckman D, Schoonhoven L, Fletcher J, Furtado K, Gunningberg L, Heyman H, et al. EPUAP classification system for pressure ulcers: European reliability study. *J Adv Nurs*. 2007;60(6):682-91.
59. Fumarola S, Allaway R, Callaghan R, Collier M, Downie F, Geraghty J, et al. Overlooked and underestimated: medical adhesive-related skin injuries. *J Wound Care*. 2020;29(Sup3c):S1-S24.
60. Blaylock B, Murray M, O'Connell K, Rex J. Tape injury in the patient with total hip replacement. *Orthopedic nursing*. 1995;14(3):25-8.
61. Brett DW. Impact on Pain Control, Epidermal Stripping, Leakage of Wound Fluid, Ease of Use, Pressure Reduction, and Cost-Effectiveness. *Journal of Wound Ostomy & Continence Nursing*. 2006;33:S15-S9.
62. Jester R, Russell L, Fell S, Williams S, Prest C. A one hospital study of the effect of wound dressings and other related factors on skin blistering following total hip and knee arthroplasty. *Journal of Orthopaedic Nursing*. 2000;4(2):71-7.
63. Koval KJ, Egol KA, Polatsch DB, Baskies MA, Homman JP, Hiebert RN. Tape blisters following hip surgery. A prospective, randomized study of two types of tape. *J Bone Joint Surg Am*. 2003;85(10):1884-7.
64. LeBlanc K, Baranoski S. Skin Tears: State of the Science: Consensus Statements for the Prevention, Prediction, Assessment, and Treatment of Skin Tears: ©. *Advances in Skin & Wound Care*. 2011;24(9):2-15.
65. Conway J, Whettam J. Adverse reactions to wound dressings. *Nurs Stand*. 2002;16(44):52-4, 6, 8 passim.
66. Widman TJO, Holly. Storrs, Frances J. Allergic contact dermatitis from medical adhesive bandages in patients who report having a reaction to medical bandages. *Dermatitis*. 2007;18(2):103.
67. Norris P, Storrs FJ. Allergic contact dermatitis to adhesive bandages. *Dermatol Clin*. 1990;8(1):147-52.
68. Lober CW, Fenske NA. Photoaging and the skin: differentiation and clinical response. *Geriatrics (Basel, Switzerland)*. 1990;45(4):36-40, 2.
69. Evans DW, Scott AT, Teasdall RD, Smith BP, Howlett A, Connor-Kerr T, et al. Mechanical properties of lower limb dermis following static and cyclic compression. *Biomed Sci Instrum*. 2012;48:104-11.
70. Parodi A, Castello M, Corbetta S, Cordera R, Cristaudo A, Cuccia R, et al. Skin and diabetes: an experts' opinion from the Italian diabetologists and dermatologists of the DiaDex group. *G Ital Dermatol Venereol*. 2018;153(5):649-58.
71. Van de Velde-Kossmann KM. Skin Examination: An Important Diagnostic Tool in Renal Failure Patients. *Blood Purif*. 2018;45(1-3):187-93.
72. Kim J, Shin Y. Medical Adhesive-Related Skin Injury Associated with Surgical Wound Dressing among Spinal Surgery Patients: A Cross-Sectional Study. *International Journal of Environmental Research and Public Health*. 2021;18(17):9150.
73. Frota OP, Pinho JdN, Ferreira-Júnior MA, Sarti ECFB, Paula FMD, Ferreira DN. Incidence and risk factors for medical adhesive-related skin injury in catheters of critically ill patients – A prospective cohort study. *medRxiv*. 2022.
74. Leu HJ. [Micromorphologic changes in the skin in primary and secondary (post-thrombotic) chronic venous insufficiency]. *Wien Med Wochenschr*. 1994;144(10-11):201-4.
75. Gao C, Yu C, Lin X, Wang H, Sheng Y. Incidence of and Risk Factors for Medical Adhesive-Related Skin Injuries Among Patients: A Cross-sectional Study. *J Wound Ostomy Continence Nurs*. 2020;47(6):576-81.
76. Bootun R. Effects of immunosuppressive therapy on wound healing. *Int Wound J*. 2013;10(1):98-104.
77. McNicol L, Lund C, Rosen T, Gray M. Medical Adhesives and Patient Safety: Consensus Statement for the Assessment, Prevention and Treatment of Adhesive-Related Skin Injuries. *Journal of Wound, Ostomy and Continence Nursing*. 2014;40(4):365-80.
78. Filimonov A, Brady JS, Govindan A, Merchant A, Eloy JA, Baredes S, et al. Postoperative complications of total laryngectomy in diabetic patients. *Laryngoscope*. 2017;127(10):2247-51.
79. Graboyes EM, Yang Z, Kallogjeri D, Diaz JA, Nussenbaum B. Patients undergoing total laryngectomy: an at-risk population for 30-day unplanned readmission. *JAMA Otolaryngol Head Neck Surg*. 2014;140(12):1157-65.
80. Nyszczy A. Total Laryngectomy Surgery in the United States and Spain: A Tale of Two Healthcare Systems: Williams Honors College, Honors Research Projects; 2019.
81. Bryant R. Types of skin damage and differential diagnosis. In: Mosby E, editor. *Acute & Chronic Wounds: Current Management Concepts*. 4th ed. 2012. p. 83-107.
82. Bryant RA. Saving the skin from tape injuries. *The American journal of nursing*. 1988;88(2):189-91.
83. Fischer H, Widdicombe JH. Mechanisms of acid and base secretion by the airway epithelium. *J Membr Biol*. 2006;211(3):139-50.

Head office:

Atos Medical AB
Hyllie Boulevard 17
SE-215 32 Malmö, Sweden
Tel: +46 (0) 415 198 00
Email: info@atosmedical.com

Manufacturer:

Atos Medical AB
Kraftgatan 8
SE-242 35 Hörby, Sweden

Always read the Instructions for use before starting to use any of the products. For Instructions for use, please visit www.atosmedical.com.

MC3093-THEN-202311

Atos
Breathing-Speaking-Living
atosmedical.com