TRACOE medical GmbH

Reichelsheimer Straße 1 / 3, 55268 Nieder-Olm / Germany

Tel.: +49 6136 9169-0, Fax: +49 6136 9169-200

info@tracoe.com, www.tracoe.com

TRACOE® care purofoam

Tracheostomy Dressing

For Adults, REF 958

For Children and Adolescents, REF 959

TRACOE® care purofoam

Trachealkompressen

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EN	Instructions for Use
DE	Gebrauchsanweisung
FR	Instructions d'utilisation
IT	Istruzioni d'uso
ES	Instrucciones de uso
PT	Instruções de uso
DA	Brugsanvisning
FI	Käyttöohje
NB	Bruksanvisning
HU	Használati utasítás
SV	Bruksanvisning
NL	Gebruiksaanwijzing
EL	Οδηγίες χρήσης
KO	사용 설명서
ZH	用于新生儿
AR	دليل استعمال

EN / Warning. Please read Instructions for Use · **DE** / Achtung. Gebrauchsanweisung beachten · **FR** / Attention. Respecter le mode d'emploi · **IT** / Attenzione. Consultare le istruzioni per l'uso · **ES** / Atención. Siga las instrucciones de uso · **PT** / Atenção. Ter atenção às instruções de utilização · **DA** / Vigtigt. Følg brugsanvisningen · **FI** / Huomio. Noudata käyttöohjetta · **NB** / OBS. Følg bruksanvisningen · **HU** / Figyelem. Olvassa el a használati útmutatót · **SV** / OBS! Se bruksanvisningen · **NL** / Let op. Rekening houden met de gebruiksaanwijzing · **EL** / Προσοχή. Ακολουθήστε τις οδηγίες χρήσης · **KO** / 주의. 사용 설명서를 준수하십시오. · **ZH** / 注意。注意使用说明书 · **AR** / تنبيه. يجب مراعاة دليل الاستعمال

EN / For single-patient use only · **DE** / Nur für einen Patienten · **FR** / Pour un seul patient · **IT** / Prodotto destinato ad un solo paziente · **ES** / Para uso por un solo paciente · **PT** / Só para um paciente · **DA** / Kun til én patient · **FI** / Vain yhden potilaan käyttöön · **NB** / Kun for én pasient · **HU** / Csak egy betegnél alkalmazható · **SV** / Endast för en patient · **NL** / Voor gebruik bij slechts één patiënt · **EL** / Μόνο για έναν ασθενή · **KO** / 단일 환자 전용 · **ZH** / 仅限于一位患者使用 · **AR** / تستعمل لمريض واحد فقط

REF **EN** / Item number · **DE** / Artikelnummer · **FR** / Référence · **IT** / Codice prodotto · **ES** / Número de artículo · **PT** / Número do artigo · **DA** / Artikelnummer · **FI** / Tuotenumero · **NB** / Artikkelnummer · **HU** / Cikkszám · **SV** / Katalognummer · **NL** / Artikelnummer · **EL** / Αριθμός είδους · **KO** / 품목 번호 · **ZH** / 产品编号 · **AR** / رقم السلعة

LOT **EN** / Batch code · **DE** / Chargencode · **FR** / Code de lot · **IT** / Numero di lotto · **ES** / Código de lote · **PT** / Código do lote · **DA** / Batchkode · **FI** / Eräkoodi · **NB** / Chargekode · **HU** / Tételkód · **SV** / Satsnummer · **NL** / Lotcode · **EL** / Κωδικός παρτίδας · **KO** / 생산 단위 · **ZH** / 批次编码 · **AR** / رمز الشحن

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For Children and Adolescents, REF 959

1. General Information

IMPORTANT: Please read the following instructions before using the product. The following instructions apply only to the above mentioned disposable dressings and must be known to the user.

Warning: The dressing should be removed if allergic reactions such as itching, redness, or other signs of hypersensitivity occur. The product must not be used in case of known sensitivity to plasters, adhesive tapes, or wetting agents.

CAUTION:

- This dressing does not contain pharmaceutical agents. Infection symptoms require medical intervention.
- The use of ointments, balms may damage or destroy parts of the dressing as well as cleaning / disinfection with hypochlorite solution or oxygenated water.

2. Product Description

The TRACOE purofoam tracheostomy dressings are made of a highly absorbent open pore polyurethane foam maintaining the appropriate moisture to wound environment. The nonsticking surface layer of the dressing side facing the skin allows quick absorption of liquids into the foam's core and features a peelable protective foil. The TRACOE purofoam tracheostomy dressing can absorb large amounts of secretions. It therefore prevents damage of the patient's skin surrounding the tracheostoma. Even in the presence of an irritated and weeping stoma the dressing does not stick to the skin, and can therefore easily be removed and replaced. The opposite side is vapor permeable and also liquid-proof when the recommended dressing change intervals are followed. This outside layer also acts as an additional mechanical, microbiological barrier. This dressing meets modern standards of moist wound treatment.

Dressing Dimensions (Width x Length x Thickness):

- REF 958 for adults: (9 x 10 x 0.5) cm, opening diameter 9 mm
- REF 959 for children and adolescents: (6.5 x 6.5 x 0.5) cm, opening diameter 5 mm.

Due to the circular opening and the zigzag slit the dressing can be easily de- or attached behind the flange of the tracheostomy tube in place. The low profile of the dressing helps to make the dressing comfortable for the patient and helps to maintain correct positioning of the tracheostomy tube within the airway.

3. Intended Use

The dressings are designed to help protect the skin surrounding and wound environment of a tracheostoma or an artificial orifice of children, adolescents and adults against abrasion (which can lead to irritation and/or infection of the skin) by:

- absorbing secretions from the tracheostoma and,
 - acting as a mechanical barrier between the skin surrounding the tracheostoma and the neck flange of the tracheostomy tube,
 - providing a waterproof but breathable (vapor permeable) barrier, that also aids in the prevention of bacterial contamination, and
 - not-sticking to the wound when the dressing is changed
- The dressing is intended for single use only.

4. Indications

The TRACOE purofoam tracheostomy dressing can be used to dress a newly established or an existing, irritated tracheostoma, as well as healed but weeping tracheostomas or artificial orifices.

5. Contra-Indications

Keep away from neonates and infants.

Do not use the dressing for infected or excessively festering wounds. Do not use for patients having diabetes mellitus or blood circulatory failure.

6. Adverse Reactions

When the identical dressing is used for too long, the dressing might cause a rash.

7. Instructions for Use

How to Use:

- Clean and dry the tracheostoma and its surrounding as usual.
- Remove the protective foil of the absorbant (skin/stoma) side. The dressing may be cut to a specific shape if necessary.
- Gently lift the tracheostomy tube neck flange to create a small gap for the dressing.
- Place the dressing with the absorbant side between skin/stoma and the tracheostomy neck flange using the zig-zag-slit. The zig-zag-slit is facing upwards.
- If necessary, the dressing can be held in place with an additional plasters or medical adhesive tapes.

How to Change the Dressing:

- Slightly lift the tracheostomy tube holding the flange and carefully pull the dressing downwards away from the patient.
- Follow again description in „How to Use“.

The dressing can absorb large amounts of liquids ($> 0.1\text{g/cm}^2$). During the early stages of treatment particularly following tracheotomy, TRACOE purofoam dressings should be inspected frequently at short intervals. Change intervals are 2 – 3 days, depending on the amount of exudation or secretions. A darker colour and swelling of the dressing, indicates that the dressing should be changed. If leakage occurs, the dressing should be changed more frequently.

8. Storage

Store in a dry location away from light and excessive humidity. Storage temperature limit is $+1^{\circ}\text{C}$ to $+30^{\circ}\text{C}$.

9. Disposal

These products may only be disposed of in accordance with the valid national regulations governing waste materials.

10. Returns

Returned used products will only be accepted subject to prior consultation and if a completed decontamination certificate is enclosed. This form is available either directly from TRACOE medical or via the website www.tracoe.com.

11. General Terms and Conditions

The sale, delivery and return of all TRACOE products shall be effected exclusively on the basis of the valid General Terms and Conditions (GTC), which are available either directly from TRACOE medical or via the website www.tracoe.com.

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