

Tracoe Purofoam-T


Comfort Meets Advanced Tracheostomy Care

A soft, breathable foam dressing with high absorbency and bacteria proof backing for advanced tracheostomy care.

0800 783 1659

info.uk@atosmedical.com

Coloplast and the Coloplast logo are trademarks of Coloplast A/S.
© 2025-11 Coloplast A/S. All rights reserved.



Tracoe Purofoam-T

The efficient way to manage tracheal secretions effectively.

Tracoe Purofoam-T dressings offer efficient, comfortable care for tracheostomy sites. Its highly absorbent foam captures secretions and exudate, while the semi-permeable film backing allows breathability, acts as a bacterial barrier, and prevents strikethrough. The soft, lightweight, non-adhesive foam minimises skin irritation during dressing changes. With a pre-cut keyhole and perforated slit, placement is quick and secure.

- ✓ Provides a rapid absorption of exudate
- ✓ Promotes moist wound healing while allowing breathability
- ✓ Helps maintain skin integrity and provides a cushioning effect around the tracheostoma
- ✓ Perforated slit design enables quick conversion from non-slit to slit dressing when needed
- ✓ Designed for extended wear - up to 7 days at the wound site
- ✓ Sterile and individually packaged



CR4555 Tracoe Purofoam-T Large



CR4557 Tracoe Purofoam-T Small

	Product name	Indication	Size	Qty
CR4555	Tracoe Purofoam-T Large	For medium to high secretion/exudate management	8 x 8 x 0.5 cm	10 pcs
CR4557	Tracoe Purofoam-T Small	For medium to high secretion/exudate management	6.5 x 5 x 0.5 cm	10 pcs

Atos (UK),
Cartwright House, Tottle Road, Nottingham, NG2 1RT
Call: **0800 783 1659**
Web: www.atosmedical.co.uk
Email: info.uk@atosmedical.com